

# Steel Frame Privacy Screen

3130/3420 Standard White

3139/3429 Sure-Chek®



ISO 9001-2008 Certified

Create a NON-CLINICAL ENVIRONMENT for your patients with Winco's Sure-Chek® colored vinyl. Style, Versatility, & Durability!

Winco's steel frame privacy screens stand up to years of use. The wide, sturdy, steel frame eliminates tipping and leaning. Both screen models come with twin wheel casters that can be converted to crutch tip bases.



## KEY FEATURES

- Wide steel frame
- Twin casters
- Replaceable screen panels

1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## STEEL FRAME PRIVACY SCREEN

3130/3420 STANDARD WHITE

3139/3429 SURE-CHEK®


- 1 Choose from our bacteriostatic vinyl panels or our OPTIONAL Sure-Chek® antimicrobial material.
- 2 Hinged panel ensures years of problem-free use.
- 3 Wide base for added stability
- 4 Smooth rolling twin wheel casters with optional crutch-tips included.
- 5 Smooth powder-coated steel frame for durability.



### SPECIFICATIONS

	3130/3139	3420/3429
Base Size	.20" x 27" (50.8 x 68.9cm)	.20 x 42" (50.8 x 106.7cm)
Overall Height	.69" (175.3cm)	.69" (175.3cm)
Overall Width	.82" (208.3cm)	.42" (106.7cm)
Single Panel Length	.49" (124.46cm)	.49.5" (125.73cm)
Single Panel Width	.27" (68.6cm)	.42" (106.7cm)
Weight	.25 lbs. (11.3kg)	.15" (6.8kg)

### COLORS



Optional

Optional

Optional

Optional

Optional

Printed colors shown may vary slightly from actual color of product.  
 Sure-Chek® vinyls are antimicrobial, antistatic, odor resistant, stain resistant, & fire retardant.