



Style, comfort, value, and great features make our Designer Care Cliner unlike any other.

The Winco 6700 premium full featured clinical recliner has fully upholstered sides that add contemporary style, modern form, and comfort for your patients. Durable and comfortably contoured urethane arms protect the upholstered finish while providing a sleek and modern style. The 6700 also features LiquiCell[®] seating that makes it more comfortable, reduces fatigue, and relieves pressure over longer periods of sitting. Optional heat and massage completes this chair model as the ultimate medical recliner.



KEY FEATURES

- Contoured urethane armrests
- Fold down side-tables
- Fully upholstered
- Trendelenburg Positioning
- LiquiCell[®] for added comfort





PRODUCT SPECIFICATIONS

DESIGNER CARE CLINER 6700

- 1 Standard Trendelenberg release on back with easy-access handle
- 2 Non-porous, fold-down side tray contains small spills. STANDARD on both sides.
- 3 Wide elastic bands on back frame and seat provide enhanced comfort and support.
- 4 **LIQUIGELL®** improves blood perfusion and aids in reducing skin pressure.
- 5 3" Nylon Total-lock and Directional-lock casters provide ease of maneuverability



SPECIFICATIONS

6700

Seat Height	.20.5" (52cm)
Seat Width	.22" (55.88cm)
Seat Depth	.19.5" (49.5cm)
Overall Height	.47" (119.4cm)
Overall Width*	.32" (81.3cm)*
Back Height Above Seat	.29.5" (74.9cm)
Overall Length (Upright)	.39" (99.1cm)
Overall Length (Reclined)	.71.5" (181.6cm)
Armrest Height	.28.25" (71.8cm)
Weight	.95 lbs. (43.1kg)
Weight Capacity	.275 lbs. (124.7kg)

*Add 10.5" to chair width for each raised side table

OPTIONS

TB-133 TB	Heat HT Heat/Massage HM	IV Pole IV	Oxygen Tank Holder O2	Versa-Table VN
Pivot-Table PN	Side-Table LT, RT, NT	Torso Support TS00	Contoured Headrest CH00	Side Cushion SC00

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



COLORS