

Padded Blood-Drawing Chair

2572/2573



ISO 9001-2008 Certified

All the features you need, and all the quality you demand.

This chair's padded seat and back offer your patients comfort. Our blood-drawing chair is available with a standard flip-up urethane arm or an optional padded vinyl pivot arm.



KEY FEATURES

- Padded seat and back
- Sturdy leveling feet
- Contoured urethane arms
- Laminated cabinet (2572 only)
- Drawer roller glides (2572 only)

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

PADDED BLOOD-DRAWING CHAIR 2572/2573

- 1 Comfortably padded seat and backrest.
- 2 Flip-up, contoured urethane armrest provides easy access and added safety for the user.
- 3 Adjustable height armrest.
- 4 High impact leveling feet ensure a sturdy chair on almost any surface.
- 5 Fully laminated cabinet and storage drawer (2572 only) with contoured urethane armrest.
- 6 Powder-coated steel frame.



SPECIFICATIONS

2572/2573

Seat Height	22.5" (57.2cm)
Seat Width	18.75" (47.6cm)
Seat Depth	18.5" (47cm)
Overall Height	44.5" (113cm)
Overall Width	2572: 35" (88.9cm)
	2573: 29" (73.7cm)
Back Height Above Seat	24.5" (62.2cm)
Overall Length	27" (68.6cm)
Adjustable Arm Height Range	28.25" - 34.25" (71.76 - 87cm)
Weight	2572: 67 lbs. (30.4kg)
	2573: 50 lbs. (22.7kg)
Weight Capacity	300 lbs. (136kg)

*Weight shown as a minimum-maximum range & varies depending on chair options ordered.

OPTIONS



TB-133 **TB**




Pivot Arm **PR**
(both models)
Dual Pivot Arm **P2**
(2573 only)



IV Pole **IV**

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



COLORS

01-26-11

REV.B

009054