

# Treatment Table 8500



ISO 9001-2000 Certified

## Customize your treatment table!

Our tables offer versatility, comfort and strength. Build onto our basic models using shelves, cabinets, and drawers to meet your specific needs.



**8500**



**8500-SH**



**8500-C1**



**8500-D1**

## KEY FEATURES

- Durable foam construction
- Sturdy steel leg frames
- Paper holder
- Adjustable leveling feet

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## TREATMENT TABLE 8500

- 1 Heavy gauge expanded vinyl & medium density foam for superior comfort.
- 2 Standard paper holder.
- 3 High strength powder-coated steel frame.
- 4 High impact leveling feet ensure a sturdy table on almost any hard surface.
- 5 OPTIONAL - Fully laminated cabinet and/or drawer



### SPECIFICATIONS

#### 8500

Overall Height*	30" (76cm) STANDARD*
Overall Width	28" (71cm)
Overall Length	72" (183cm)
Weight**	74 lbs. (33.6kg)**
Weight Capacity	400 lbs. (181kg)

\*Special height legs available in 19", 25", & 36"

\*\*Weight varies depending on table options ordered.  
Add 25 lbs. for 8500-SH

### OPTIONS



TB-133 TB



\*Special Height Legs




Shelf SH



Cabinet **C1**  
Double Cabinet **CC**  
Cabinet & Drawer **CD**




Drawer **D1**  
Double Drawer **DD**




Matching Pillow **MPOO**

#### STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

#### OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



#### SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



### COLORS