

# Adjustable Back Treatment Table

8570



ISO 9001-2000 Certified

Superior durability, outstanding features, and an affordable price

Our adjustable back tables offer the durability and range of standard features you expect at a price that you'll love. One look, and you'll see the difference.



**8570**

## KEY FEATURES

- Gas spring adjustable backrest
- Durable foam construction
- Sturdy steel leg frames
- Paper holder
- Adjustable leveling feet



**8570-SH**



**8570-C1**



**8570-D1**

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## ADJUSTABLE BACK TREATMENT TABLE 8570

- 1 Adjustable headrest to enhance patient comfort. Angle range: 6-57 degrees
- 2 One-touch handle easily raises headrest to desired angle.
- 3 Standard paper holder.
- 4 High strength powder-coated steel frame.
- 5 High impact leveling feet ensure a sturdy table on almost any hard surface.
- 6 Heavy gauge expanded vinyl & medium density foam for superior comfort.



### SPECIFICATIONS

#### 8570

Overall Height*	30" (76cm) STANDARD*
Overall Width	28" (71cm)
Overall Length	72" (183cm)
Adjustable Backrest Range	6 - 57 degrees
Weight**	91 lbs. (41.3kg)**
Weight Capacity	400 lbs. (181.4kg)

\*Special height legs available in 25" or 36"  
 \*\*Weight varies depending on table options ordered.  
 Add 25 lbs. for 8570-SH

### OPTIONS



TB-133 TB



\*Special Height Legs



Shelf SH



Cabinet C1




Drawer D1




Matching Pillow MPOO

### STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

### OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



### SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



### COLORS